

FIXATION & LACZ STAINING OF MOUSE EMBRYOS (BLACK LAB)

1. Euthanize female mouse
 2. Remove uterus to a petri dish of PBS-A
 3. Remove embryos with yolk sac
 4. Remove all membranes (**including the amnion**); save yolk sac in tail buffer for Southern, if necessary
 5. Rinse embryos in PBS-A to remove all blood, etc.
 6. Fix embryos in FIX SOLUTION at 4°C in 24- or 6-well plates (5 ml fix in a 6 well; 1 ml fix for 24 well)

8.5-9.5 dpc	30 min.
10.5-11.5 dpc	60 min.
12.5 dpc	1.5 hours
13.5 dpc	2 hours
14.5 dpc	2.5 hours
15.5 dpc	3 hours
16.5-17.5 dpc	4 hours

For 15.5 dpc and younger, no skinning is required; for 15.5 and older, skinning is required
 7. Rinse quickly with PBS-A; then incubate in PBS-A for 30 min.
 8. Put embryos in STAINING SOLUTION overnight in the dark (wrap in foil) at RT
For 13.5-15.5 dpc embryos, the 'bros should be permeablized with detergent during staining
make staining solution FRESH
 9. Rinse quickly with PBS-A; then incubate in PBS-A for 30 min.
 10. FIX again overnight in 4% formaldehyde (salmon colored--see below)
 11. Embryos can be stored indefinitely in 4% formaldehyde or in PBS-A at 4°C
-

SOME EMBRYOS MAY NEED TO BE CLEARED FOR WHOLE MOUNT PHOTOGRAPHY

(especially true for 12.5 and older; skinned embryos CANNOT be cleared)

1. following fixation, rinse embryos in PBS-A
2. DEHYDRATE embryos in ethanol
 - put embryos in GLASS scintillation vials
 - add 10 ml 50% ethanol; incubate for 30 minutes
 - add 10 ml 70% ethanol; incubate for 30 minutes
 - add 10 ml 100% ethanol; incubate for 30 minutes
 - add 10 ml 100% ethanol; incubate OVERNIGHT (can be stored at -20°C indefinitely)
3. Next day, add 10 ml 100% ethanol; incubate for 30 minutes
4. Add 10 ml 1:1 benzyl benzoate: benzyl alcohol; WATCH CAREFULLY, embryos will clear quickly
5. To photograph put in glass petri dishes in FRESH 1:1 benzyl benzoate: benzyl alcohol
be sure that embryos are completely submerged for photography

ALL WORK FOR CLEARING MUST BE IN **GLASS** AND MUST BE **ANHYDROUS**

(plastic will melt and water will cause a white precipitate that cannot be removed)

SOLUTIONS

10% formaldehyde from paraformaldehyde

mix 50 ml dH₂O with 10 g paraformaldehyde

heat and stir on a hot plate (solution will remain cloudy)

add ~100 µl of 1.0 N NaOH which will make the solution go clear

filter the solution through a 0.8 micron filter (0.45 micron filter can also work)

add 10 ml 10X PBS-A

QS to 100 ml with dH₂O

aliquot and store at -20°C

FIX SOLUTION for 10 ml

2% formaldehyde 2 ml of 10% stock

0.2% glutaraldehyde 80 µl 25% stock (Sigma EM grade)

1X PBS-A 1 ml of 10X stock

6.92 ml of dH₂O

add 5 µl of phenol red (solution should be SALMON colored--light pink)

if necessary, adjust pH by adding 10% HCl until proper color (a few microliters)

STAIN SOLUTION (10 ml)

0.5 ml 100 mM potassium ferrocyanide (store stock in 1 ml aliquots at -20°C)

0.5 ml 100 mM potassium ferricyanide (store stock in 1 ml aliquots at -20°C)

20 µl 1.0 M MgCl₂

100 µl Xgal in DMF (100 mg/ml) (store stock at -20°C)

1.0 ml 10X PBS-A

7.9 ml dH₂O

To permeablize--

add NP40 to a final conc of 0.02%

20 µl of 10% NP40

add sodium deoxycholate at 0.01%

20 µl of 5% sodium deoxycholate

add Tris pH=7.4 to a conc of 20mM

200 µl of 1.0 M Tris pH=7.4

10X STOCK OF PBS-A (500 ml)

NaCl 40.0 g

KCl 1.0 g

KH₂PO₄ 1.0 g

Na₂HPO₄ 5.75 g

dissolve into a final volume of 500 ml dH₂O and autoclave